Nyano Nani न्यानो नानी

User Manual

Version 2.0


Table of Contents

Safety Guidelines	1
Technical Specifications	5
Theory of Operation	6
Operation Instructions	7
Warning, Operation & Troubleshooting	10
Preventive Maintenance	14

The user should carefully read these safety and usage instructions before the use of the device.

This device is not suggested to be used in any areas where combustible and explosive gas mixtures are present.

Electrical connections to equipment which is not listed in these instructions for use should only be made following consultation with the respective manufacturers.

WARNING

Disclaimer

The design of the equipment, the accompanying literature, and the labeling on the equipment take into consideration that the purchase and use of the equipment are restricted to trained professionals, and that certain inherent characteristics of the equipment are known to the trained operator.

This publication excludes references to various hazards which are obvious to a medical professional and operator of this equipment, to the consequences of product misuse, and to potentially adverse effects in patients with abnormal conditions.

Maintenance

The device must be inspected and serviced regularly by trained and authorized personnel every 4 months. Any necessary repair and replacement of spare parts should be done by authorized personnel only. Failure to do so may compromise the correct functioning of the device.

Caution

It is the responsibility of the appropriate clinical decision-maker to draw conclusions from the skin temperature measurements taken.

Accessories

Do not use any spare accessories other than those provided by authorized service personnel. Product modification or misuse can be dangerous.

Safety Instructions

- It is advisable to not leave an INFANT unattended under the INFANT WARMER.
- Radiant warming increases insensible water loss. Appropriate measures to maintain proper fluid balance should be considered.
- Permissible distance between the radiant heater and the bottom of the bassinet is 75cm. No changes to this distance should be made by any means since a larger distance between the two points may cause insufficient heat on the INFANT and a reduced distance may result in excessive heating.

Safety Instructions (Contd...)

- The skin temperature sensing probe should be in direct contact with the skin to provide accurate monitoring of the infant's skin temperature. Failure to maintain direct skin contact can result in overheating and possible burning. Check the INFANT's condition at least every 15 minutes for correct Sensor attachment, reddened skin areas, and proper skin temperature.

Note: Rectal temperatures are not appropriate for controlling the heater output.

- Do not leave the infant unattended in the Cart/Bassinet when the side panels on the front panel are folded down.
- Inspect the side panels of the bassinet to prevent INFANT falling out.
- The bassinet withstands a maximum of 20 kg load. It is recommended to operate the INFANT warmer with 50% of the maximum load at all times.
- External heated mattress, use of phototherapy or sunlight can cause an increase in the temperature of INFANT abruptly.
- The INFANT WARMER should be positioned on its central axis and the bassinet should not be tilted to horizontal or vertical sides. This may affect the performance of the INFANT WARMER.
- INFANT WARMER is not suitable for use in presence of flammable anesthetic gases or other flammable materials, such as some cleaning fluids.
- Environmental conditions like AIR movement, Humidity can affect the thermal (temperature) balance of the INFANT and may take some time to be stable.
- The INFANT WARMER shall be used only by appropriately trained health personnel and under the direction of qualified medical personnel who are familiar with currently known RISKS and BENEFITS of warmer use.

Definitions

WARNING: A WARNING statement contains important information about a potentially hazardous condition which can result in death or a serious injury to the patient if not avoided.

CAUTION: A CAUTION statement contains important information about a potentially hazardous condition which can result in a minor or even severe injury to the patient if not avoided.

DEVICE: An infant radiant warmer with the name “Nyano Nani” designed and developed completely in Nepal.

NIC: National Innovation Center, Nepal. The manufacturer for this DEVICE.


LIMITATION OF LIABILITY

ALTHOUGH THE DEVICE HAS BEEN DESIGNED AND TESTED WITH UTMOST PRIORITY GIVEN TO THE SAFETY OF THE PATIENT AND TO MINIMIZE ANY HAZARDS, THE MANUFACTURERS (NIC) OR THE FUTURE COMPANY ASSOCIATED WITH THE DEVICE, IN NO EVENT SHALL BE LIABLE FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, PUNITIVE OR CONSEQUENTIAL DAMAGES TO THE PATIENT AND/OR THE DEVICE WHATSOEVER RESULTING FROM ANY (I) USE OF THE DEVICE, (II) FAILURE TO FOLLOW THE SAFETY AND USAGE INSTRUCTIONS, (III) INTERNAL ACCESS, MODIFICATION OR TAMPERING TO THE DEVICE, (IV) ERROR AND ANY INACCURACIES ON THE CONTENT PROVIDED BY THE MANUFACTURERS INCLUSIVE BUT NOT LIMITED TO THIS DOCUMENT.

Technical Specifications

Power Requirements	220-240V AC System 50 Hz 700 Watt.
Overload Protection	5A Circuit Breaker
Examination Light	3 Watt. LED
Heater	Elstein Ceramic 650 Watt.
Sensors	Skin & Air PT100
Coating	MS Powder Coating
Reflectors	Parabolic Steel Reflectors
Temperature Controller	Auto-PID
Bassinet	Acrylic Bed Adjustable Side Panels
Electrical Safety	Earthing with Insulated Body

Theory of Operation

•••

The radiant heat is directed by parabolic reflector onto the infant. The warmer is controlled by a single heater driver with automatic skin temperature control.

Examination (Observation) light in the warmer module provides added illumination of mattress area.

Operation Instructions

•••

1

Place the baby warmer on a levelled surface and make sure **the wheels are locked**.

2


Plug In the power supply to the wall and **Turn On** the switch on the side of the device.

3

Warmer will automatically go into **Pre-Heat** state, as indicated by the display.

Warming Time: ~20 mins.

4

Place the baby in the **Bassinet** after the **Buzzer** and attach the **Skin Sensor** on a preferred location.

5

Press **SET** to program desired temperature and use the **+** and **-** buttons to adjust **SET TEMPERATURE**.

Press **SET** again to exit.


Warnings, Alarms & Troubleshooting

Note: Alarm rings for power failure, high skin and air temperature, low skin temperature and sensor failure.


Power Fail

Cause	AC Power to the unit is interrupted.
Indication	Display shows the Pr FAIL message with continuous buzzer.
Measures	Check the AC power supply. After power is restored, the alarm will stop automatically.

Note: The sensors are attached to the front panel with fasteners, please ensure that the connector is tightly secured.

If the **skin sensor fails**, indicator next to the skin temperature display will **glow** and vice versa for the air sensor.


Sensor Fail


Cause	Temperature sensors are unplugged
Indication	Display shows the SEnS FAIL message with continuous buzzer.
Measures	Check if the sensor connectors are loose. After they are adjusted, the system will resume its operation.

Note: Air temperature sensor must be placed on the bassinet in open air, free from obstructions, facing upwards. Please make sure that the metallic face of the sensor is not covered by the blanket or has not come into contact with the metal & acrylic frame of the bassinet.

High Air Temperature

Cause	Air temperature is higher than 38.5°C .
Indication	Display shows the HIGH Air message with continuous buzzer.
Measures	Make sure the air sensor is not attached to the baby. Heater power indicator should be OFF .

Low Skin Temperature

Cause Skin temperature is lower than set point.

Indication Display shows the **LO tP** message with continuous buzzer.

Measures Check for the heater indicator, it should be **ON**.

High Skin Temperature

Cause Skin temperature is higher than set point.

Indication Display shows the **HIGH tP** message with continuous buzzer.

Measures Check for the heater indicator, it should be **OFF**.

Preventive Maintenance

✓ Cleaning the Equipment

The equipment must be thoroughly cleaned and disinfected **once every week** and **every time the infant is removed from the bassinet**.

Before and after each use, clean the **sensor tip** with **alcohol swab**.

✓ Sensor Maintenance

Ensure the sensors are **properly fastened** in the **sensor slots**.

Do not **twist or jerk** the sensor wires.

✓ Heater Maintenance

Remove the heater hood from the top section and **check the heater wires** **routinely**.

Preferably **every two months**.

Scan this code to **download**
the **digital version**

Note: The product satisfies all the guidelines provided by IEC 6061-2-21 to ensure maximum safety to the child.

**Designed & Developed @ National
Innovation Center Nepal**

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